



**Med  
Page**

Celebrating  
35 Years of  
Innovation

Quality Care Solutions  
**Easylink UK**

## CARING SOLUTIONS 2019/2020

*Assistive technologies to aid independent living*



FM 84603



3 Melbourne House | Corby Gate Business Park | Corby | Northants | NN17 5JG  
T +44 (0) 1536 264 869 | F +44 (0) 1536 269719

**[www.medpage-ltd.com](http://www.medpage-ltd.com) | [www.easylinkuk.co.uk](http://www.easylinkuk.co.uk)**

[f http://facebook.com/MedPageLtd](https://www.facebook.com/MedPageLtd) [@MedPageLtd](https://twitter.com/MedPageLtd)

Medpage Limited T/A Easylink UK. Company registered in England & Wales No. 05468599

## Visit our websites to view our entire product range

<b>Pages 3 - 5</b>	Seizure movement detection alarms
<b>Page 6</b>	Observation cameras (RF & Wi-Fi)
<b>Page 7</b>	Brief advice on product selection
<b>Page 8</b>	Memory & cognitive impairment aids
<b>Page 9</b>	GPS location trackers
<b>Page 10</b>	Daily living aids
<b>Page 11</b>	Medpage long-range transmitter products
<b>Page 12</b>	Medpage MPPL range alarm receivers
<b>Page 13</b>	Medpage MPPL range transmitters
<b>Page 14</b>	Rondish brand alarm receivers
<b>Page 15</b>	Rondish brand transmitter products
<b>Pages 16 - 21</b>	Call pendant and special switch products
<b>Pages 22 - 24</b>	Door alarms for wander detection
<b>Pages 25 - 27</b>	Movement sensing alarms (PIR's)
<b>Pages 28 - 29</b>	Floor pressure mat alarms
<b>Pages 30 - 32</b>	Bed leaving detection alarms
<b>Pages 33 - 35</b>	Chair leaving detection alarms
<b>Page 36</b>	Off-site alert home care alarm systems

E&OE. We endeavour to maintain 100% accuracy with the content of our publications. We reserve the right to alter pricing and product specification without notice. Please check our websites for any changes.

**Trademarks:** Medpage and MemRabel are Trademarks owned by Medpage Limited.

**Orders:** We accept orders from the public, Local Authorities, NHS and trade resellers (subject to account approval).

**Terms and Conditions - Returns Policy - 28 Day Product Trials - Warranties, GDPR your personal data privacy** Details available on our website.

# EPILEPSY - MEDPAGE SEIZURE MOVEMENT DETECTION ALARMS

## MP5 Seizure movement detection monitor with pager alarm

The MP5 was first introduced to market in 1995 and has remained our top selling Tonic/Clonic seizure monitor ever since. A patented movement sensor is positioned under the person's mattress to detect nocturnal Tonic/Clonic seizure movement. An alert is transmitted to carer pagers (2 supplied) upon detected seizure movement.

- Wireless alerts signalled to MPPL pager with tone or vibration alert
- Adjustable sensitivity to suit patient body weight and bed mattress type
- Adjustable movement delay alarm virtually eliminates false alarms
- Digitally stored settings prevent tampering
- Mains powered with 36-hour battery backup
- Built in microphone detects vocalisations preceding or during a seizure
- Supplied with 2 carer pagers (1 free)
- Can be integrated into existing Nurse Call and Telecare systems



**ORDER CODE MP5 PRICE £205.00 (£246.00)** VAT Exemption available

## MP5-UT Children's epilepsy seizure movement detection alarm - Ultra-Sensitive

The MP5-UT uses patented sensor technology specifically engineered to detect Tonic/Clonic seizure movement while a child is in bed.

- Detected seizures are notified via the MPPL pagers with tone or vibration alert
- Suitable for babies from age 1 month to adults
- Supplied with 2 carer pagers (1 free) so that both parents are alerted to detected seizures
- Digital sensitivity control
- Sound sensor detects vocalisations preceding or during a seizure
- Digital alarm delay control virtually eliminates false alarms by ignoring natural sleeping movements such as turning over
- Tamper proof settings & fault reporting to pagers



**ORDER CODE MP5-UT PRICE £170.00 (£204.00)** VAT Exemption available



## MP5-UTB Epilepsy bed movement seizure detection alarm

The MP5-UTB is a reliable nocturnal seizure movement detection alarm for babies through to adults at an affordable price. A portable battery-operated alarm with flashing coloured lights and adjustable audible alarm raises an alert to carer of detected seizure movement.



- Adjustable sensitivity to suit the needs of the user
- Digital alarm delay control virtually eliminates false alarms by ignoring natural sleeping movements such as turning over
- Digitally stored tamper proof settings
- Inbuilt sound sensor detects vocalisations preceding or during a seizure
- Optional plug in under pillow vibrating pad to assist deaf carers (VIB-PAD)

**ORDER CODE MP5-UTB PRICE £170.00 (£204.00)** VAT exemption available  
See page 12 for optional power supply and vibrating pillow shaker for alarm unit

## MP2-V3K Medpage seizure movement detection alarm for care homes

The Medpage MP2-V3K system is designed to offer multiple patient seizure movement detection. The MP2-V3 monitor uses advanced sensor and computer software processing to positively identify and raise alarms to ongoing seizures. The PAG-11c pager can be used solely as the alarm device or in conjunction with the CMD-11 Desktop/Wall mounting call station. Active seizure alarms are clearly displayed on the pager and call station with user identity, and sensor type if required.

- Wireless alarm to carer message pager and/or CMD-11 & ILB-11 (see P14)
- Suitable for use in larger buildings
- Adjustable sensor sensitivity
- False positive alarm delay control
- Suitable for babies from 1 month to adults
- Detection of Tonic/Clonic seizure movement
- Vocalisation microphone with alarm to pager
- Can connect to Nurse Call and/or Telecare equipment



Seizure movement monitor system including 1 x PAG-11c pager

**ORDER CODE MP2-V3K PRICE £212.00 (£254.40)** See website for more details

# EPILEPSY - MEDPAGE SEIZURE MOVEMENT DETECTION ALARMS

## MP5-ULTRA Advanced seizure movement detection alarm with detection mode selection and carer alarm pagers

The most advanced epilepsy monitor in the MP5 range. Advanced software algorithms continually analyse sleeping movements to positively identify seizure activity. Suitable for complex epilepsy and specialist beds including airflow mattresses.



- Detected seizures are notified via MPPL pager (2 supplied) with tone or vibration alert
- Specialist customisation available
- Programmed with selectable seizure detection modes
- Digital movement delayed alarm control
- Inbuilt sound sensor to detect vocalisations preceding or during a seizure
- Large area patented monitoring sensor technology offers improved movement detection
- Digital alarm delay control virtually eliminates false alarms by ignoring natural sleeping movements such as turning over
- Tamper proof settings & fail safe alarm to pagers
- Tonic/Clonic and seizure movement with twitching
- Myoclonic cluster seizures
- Airflow mattress mode
- All Medpage seizure movement monitors have 36-hour battery backup
- Suitable for all ages: babies - adults
- Nurse Call and Telecare systems output connection

**ORDER CODE MP5-ULTRA PRICE £439.00 (£526.80)** VAT exemption available

**All Medpage seizure movement monitors are supplied with a 36 month warranty**

Medpage Limited T/A Easylink UK, established 1984, have provided solutions to aid epilepsy care for over 25 years, specialising in the design and manufacture of monitors that detect seizure movements from a sleeping person in bed. To ensure the highest product quality and safety, our company is accredited by BSI to ISO 9001:2015.

At Medpage Limited, we take great care in the design, manufacture and after sales support for our seizure movement alarms. When you purchase an epilepsy seizure movement alarm from Medpage, our highly trained team are readily available for help, advice, and technical support, if needed, to ensure your monitor operates at maximum efficiency. Our technical support is provided at no extra cost.

## OBSERVATION CAMERAS (RF & Wi-Fi)

### KFM30C Recording Wi-Fi Camera for home and business

View live video monitoring direct to your Smartphone with this simple to install and use home or business WiFi connected camera.

- Audio & video recorded to SD card (sold separately)
- 360° camera rotation via Smartphone APP
- Hi-Def night vision IR LED sensors
- Link with trigger sensors - PIR, Door, etc.
- Powered by 5V AC/DC adaptor (included)



Complete camera system **ORDER CODE KFM30C PRICE £39.99 (£47.99)**

### MEDBM-05 Wireless digital baby monitoring camera system

Useful camera system ideal for baby monitoring and patient care.

- Portable rechargeable video monitor
- Battery standby 4-6 hours (typical)
- DC adaptor powered camera
- Wireless range up to 60M (indoors)
- Multiple camera (up to 4) option
- Camera Pan/Tilt from monitor
- High performance night vision



Complete camera system **ORDER CODE MEDBM-05 PRICE £83.32 (£98.98)**

### MEDBM-03 Digital video watch wireless baby monitor

Hi performance digital wireless video monitor camera with watch style live video viewer.

- VOX vibration alert for camera detected sounds (ideal for deaf carers)
- Wireless range up to 60M indoors (typical)
- Camera night vision viewing
- 2-way voice intercom



Complete camera system **ORDER CODE MEDBM-03 PRICE £94.00 (£112.80)**



## A proud history of innovation and passionate service

Medpage Limited has manufactured senior and specialist care products for over 35 years and is accredited by BSI to ISO 9001:2015.

Many of our products are developed from personal family care scenarios requiring Assistive Technology, even more so recently with prevailing falls and Dementia.

The majority of our products use radio transmission technology to relay signals from a sensor to an alarm device. All radio signalling products achieve RED compliance and are certificated where required.

Quality and performance at an affordable price has remained our mission for 35 years. You can rest assured that you will receive high quality service and courtesy from our highly trained team.

In the following pages you will find a variety of products, which can help with:

**Memory prompting** - Automated reminders for medications, chores etc.

**Reducing the risk of falls** - Worn fall detectors, bed and chair occupancy alarms

**Wandering away from home** - GPS location products

**Wandering in the home** - Movement sensors and door/window alarms

**Personal help call buttons** - Alert a carer wirelessly that help is needed

You will find various options for achieving an alarm solution. This may be the type of transmitter or type of alarm device. Some factors to consider are:

### **Domestic, Care Home or Hospital use**

A single or multiple application (example: 1 door alarm or alarms for multiple doors)

Transmitter range (based on size or construction of building)

**For single alarms** - The MPPL pager or NMDRX tabletop alarm is a popular choice

**For up to 5 different transmitter devices** - The NMDRX alarm station is suitable

**For mixed long range alarms/sensors** - The POCSAG digital message pager would be the best option

**For multiple-patient and mixed device monitoring for Care Home or Hospital** -

Rondish systems provide more flexibility and integration for larger buildings

**For off-site warning alerts** - You may consider an automated telephone dialler or Wi-Fi connected product

Our website contains the full product specification for all items detailed in this brochure. Simply visit our website and enter the product code into the search box.

If you need any assistance our technical advisors are available to offer advice and technical support if required.

# MEMORY PROMPTING & COGNITIVE IMPAIRMENT AIDS

## MEMRABEL 2 Digital memory prompting calendar clock

This unique clock allows a family member or health visitor to set daily timed memory prompts using voice, picture or video to plan a person's daily routine.

- Create personalised memory prompts on your phone or PC, or choose from over 100 pre-installed reminder videos and graphic/text prompts
- File transfer via flash drive or SD card
- Select reminder playback for one-off, daily, weekly
- Choice of time display including traditional analogue
- Auto seasonal time adjustment
- Powered by mains AC/DC adaptor



**ORDER CODE MEMRABEL2ND PRICE £99.99 (£119.00)**

## MEMRABEL 2i Create memory prompts from your Smartphone

All the features of the Memrabel 2 but, you can now use your Smartphone to set up and send video, audio or picture memory prompts with text that will display on the MemRabel 2i at times you choose.

- Program reminders from anywhere via your Smartphone
- Aid for Stroke, cognitive impairment, learning disability, dementia, Autism
- User friendly Smartphone APP iOS and Android
- Requires home broadband
- 48 Months free subscription £58 P/A thereafter



**ORDER CODE MR2 PRICE £166.66 (£199.99)**

## MED-01V Talking medication reminder alarm clock

Versatile talking reminder alarm clock that speaks the time and date at the press of a button, is also a handy daily memory prompting aid. Record up to 6 personalised alarm messages.

- Medication reminders with daily repeat
- Appointments
- Daily chores, eating and drinking reminders
- Powered by 2 x AA batteries or optional AC adaptor



**ORDER CODE MED-01V  
PRICE £29.95 (£39.99)**



## GPS LOCATION TRACKERS - QUICKLY FIND A MISSING PERSON

### MMTBO2 Medpage Micro-GPS tracker with fall sensor & secure Smartphone APP

Rest assured knowing that should a person go missing, they can be located immediately with the Medpage micro tracker.

- GPS location accurate to 5M
- Geo-fence safe zone with breach alarm via APP
- Inbuilt fall sensor sends an alert to emergency contacts should the user experience a fall
- Location by GPS/Wi-Fi hub/GPRS
- Two-way mobile voice calling
- No monitoring or subscription charges
- SOS button calls carer and sends location details
- Includes lanyard and wrist strap



Requires mobile SIM **ORDER CODE MMTBO2 PRICE £74.99 (£89.99)**

Mobile Roaming SIM offer. Telephone sales for details

### MED-V36 Medpage GPS waterproof tracker watch with heart rate & BP monitor

Lightweight GPS personal tracker watch offers discrete peace of mind to those caring for a loved one who is prone to wandering.

- GPS tracking accurate to within 5M
- Send personalised reminder alarms to the watch from carer Smartphone - watch vibrates with tone alarm to alert user
- Approximate Heart Rate and Blood Pressure check from carer Smartphone via control APP
- Geo-fence & Wi-Fi zone protection with breach alarm
- SOS function sends location text & 2-way voice calling
- No subscription or monitoring charges
- Waterproof to 1M (for 30 minutes maximum)
- Requires mobile SIM (Optional roaming SIM deal)



**ORDER CODE MED-V36 PRICE £78.20 (£93.90)**

### MemRabel Remote controlled wireless object locator COL-02

People in the early stages of dementia frequently misplace items around the home. Now, it's very easy to find them. Place a remote beeper on the object (6 included) then simply press the corresponding button on the remote to activate the colour coded beeper.

- Includes 6 beepers powered by CR2032 battery
- Attach to keys, handbags, medicines etc.
- Suitable for the visually impaired - 50M location range



**ORDER CODE COL-02 PRICE £18.29 (£21.95)**

### GMRC Easy to use TV remote control - simplifies control of TV

Most of the remote controls supplied with TV's and audio players are extremely complicated in appearance and use. This compact remote makes it very easy for anyone, especially those who may be confused.

- Learns control codes from existing remote controls
- Can control 2 devices (TV & Set box) using buttons A or B
- Powered by 2 x AAA batteries (included)
- Can be used with existing remotes
- Useful to aid people with poor dexterity
- Ideal for people with sight impairment or misplaced glasses



**ORDER CODE GMRC PRICE £10.82 (£12.98)**

### MAMR-02 Motion activated recordable message player

Simply record a personalised reminder (up to 12 seconds), switch on the reminder player, when a person walks into the detection beam the recording is played back.

- Useful for dementia and cognitive impairment care
- Supplied with wall fixing swivel bracket
- Powered by 3 x AA batteries (included)
- Beam detection range 2-3M
- 3 minute delay between movement triggers



**ORDER CODE MAMR-02 PRICE £16.62 (£19.94)**

### YBZP PIR motion activated battery powered safety light

Useful elderly care aid. Provides a good source of light for bed leaving etc.

- Home safety lighting - handy PIR battery light
- 8 LED PIR sensor light
- Powered by 3 x AA batteries
- Direction adjustable lighting (up/down)
- Fixings include wall bracket and Velcro dots
- On/off switch with manual and auto-sensor mode
- Dims: L 160mm W 50mm D 20mm



**ORDER CODE YBZP PRICE £8.29 (£9.95)**

## Long-Range POCSAG data message encoded transmitter products

The products detailed on this page are featured in various kits in the following pages. An unlimited number of transmitters can be used with single or multiple pagers. The messages are factory programmed so you need to provide us with your messages when placing an order for additional or multiple products. Single items are shipped with the appropriate default message. Some items are not featured in this brochure.



**POC-840DC** Quick and easy to install door & window alarm with on/off switch. Programmable with message up to 40 characters.

**ORDER CODE POC-840DC PRICE £50.00 (£60.00)**

400M range  
Power 1 x CR2032  
Low battery warning



**POC-940** Waterproof (IP67) call pendant with wrist strap and lanyard included. Programmable with message up to 40 characters.

**ORDER CODE POC-940 PRICE £46.00 (£55.20)**

400M range  
Power 1 x CR2450  
Transmit confirmation LED - Low batt alert



**NMDTX** Multi-function message programmable POCSAG transmitter with configurable N/O N/C input.

**ORDER CODE NMDTX PRICE £51.00 (£61.20)**

400M range  
Power 2 x AAA  
Bed/chair & floor pressure mat modes.



**EM300** Digital data message pager with tone/vibrate alert. Stores last 20 pages with time & date.

**ORDER CODE EM300 PRICE £61.00 (73.20)**

Sturdy belt clip  
Power 1 x AAA  
Time/date calendar  
Time alarms



**EM400** Digital data message pager with powerful vibration and tone alert. Top view vertical display.

**ORDER CODE EM400 PRICE £56.00 (£67.20)**

400M range  
Power 1 x AA  
LCD backlit display  
Message review



**BC346** Long range call button with inertia and gravity fall sensor (software adjustable) IP67 waterproof.

**ORDER CODE BC346 PRICE £78.00 (£93.60)**

400M range  
Power 1 x CR2450  
Message programming  
Transmit & fall LED



**EM1600U** Signal repeater for all Medpage long-range transmitter products (this page). Extend range up to 500M. Suitable for Care Home and security applications

**ORDER CODE EM1600U PRICE £135.00 (£162.00)**

Range up to 500M  
Power 12V DC @ 1A



**MP-900** Radio frequency to Wi-Fi adaptor. Receive alerts to pagers and Smartphone Via free APP

**ORDER CODE MP-900 PRICE £95.00 (£114.00)**

Requires home or business Broadband.  
Perfect combination for security and home care.



## MEDPAGE WIRELESS ALARM RECEIVERS & ACCESSORIES

Medpage manufacture the alarm receivers detailed below to offer choice for alerting carers. Some people prefer a pager, others a tabletop alarm. The CTMV allows recording of personal alarms or reminders, triggered by remote transmitter. Where you see the Medpage logo in the following pages, it indicates the products are compatible with these alarm receivers.



Portable alarm pager with belt clip. Robust in design and very easy to use. Compatible with all MPPL transmitters.

Power: 2 x AA  
Tone/vibration alarm  
Adjustable volume  
For up to 9 devices  
Typical 100M range

**ORDER CODE MPPL2**  
**PRICE £31.00 (£37.20)**



Portable 5-Channel alarm station with tone & individually coloured flashing light alarms. Compatible with all MPPL transmitters.

Power: 3 x AA  
Tone/flashing alarm  
Adjustable volume  
For up to 5 devices  
Typical 100M range

**ORDER CODE NMDRX**  
**PRICE £37.49 (£44.99)**



Portable alarm receiver with 3 recordable voice channel alarms of up to 15 seconds. Relay output allows any Medpage transmitter to connect to Nurse Call or Telecare alarm.

Power: 3 x AAA  
Voice recorder alarm  
Adjustable volume  
For up to 3 devices  
Typical 100M range

**ORDER CODE CTMV**  
**PRICE £35.83 (£42.99)**



Plug-in vibrating under-pillow shaker for use with the NMDRX alarm station for the deaf and hard of hearing.

Rated at 300mA  
Cord length 2M  
3 Volt (MAX 5V)  
3.5mm Jack plug

**ORDER CODE VIB-PAD**  
**PRICE £12.00 (£14.40)**



Medpage signal range extender compatible with all MPPL/Medpage transmitters.

Power: 5V DC adaptor  
36 hour battery backup  
Typical 100M range  
Auto transmitter pairing

**ORDER CODE MPPL-RPT**  
**PRICE £84.00 (£100.80)**



UK Plug top AC/DC power adaptor  
Please enquire for correct product model

## MEDPAGE BRAND WIRELESS TRANSMITTER ALARMS

All products featured in this section are additional transmitters for use with the MPPL pager, NMDRX 5-Channel alarm or CTMV 3-Channel alarm. They can be ordered separately where extra or mixed sensors are required. All products have an approximate transmitter range of 100M on 433MHz. Use the product code to search our website.

Universal transmitter with N/O N/C volt free trigger input

**ORDER CODE NMDTXM PRICE £33.32 (£39.99)**



Waterproof push button call pendant with lanyard and wrist strap

**ORDER CODE MPPL-ERTX Price £14.00 (£16.80)**



Wireless door/window alarm transmitter

**ORDER CODE MED-DCT PRICE £11.66 (£13.99)**



Wireless doorbell button transmitter

**ORDER CODE NMD-DB1 PRICE £8.28 (£9.94)**



Sound activated transmitter for baby monitoring and voice help call

**ORDER CODE MPPL-SAT Price £21.00 (£25.20)**



Battery operated PIR movement sensor transmitter

**ORDER CODE MED-PIR2 PRICE £14.16 (£16.99)**



Wrist worn fall sensor with touch sensor call button

**ORDER CODE MPPL-FA6 PRICE £27.45 (£32.94)**



Washable cotton enuresis sensor with alarm transmitter - Long range version available. Please enquire for details.

**ORDER CODE NMDTXM-ES01 PRICE £61.72 (£74.06)**



PIR Movement sensor transmitter with external ON/OFF switch

**ORDER CODE MPPL-MS PRICE £28.00 (£33.60)**



# RONDISH BRAND WIRELESS ALARM RECEIVERS



Medpage is the UK distributor for Rondish nursing care products. The products featured on this page are alarm receivers, which can be used individually or in combinations for individual or multiple room/ward monitoring care systems. Where you see the Rondish logo in the following pages, it indicates the products are compatible with these alarm receivers.

## ILB-11 Wall mounting audible & visual alarm



- Can pair with up to 32 transmitter products
- Audible alarm with volume adjustment
- Bright flashing alarm LED's
- Integral signal booster for all paired transmitters

**ORDER CODE ILB-11 PRICE £112.00 (£134.40)**

## CMD-11 Desktop alarm station with signal booster



- Desktop wireless caller display alarm (99 devices)
- Integral signal booster for paired devices
- Mains powered with battery backup (10 hours)
- Multiple call display - call review

**ORDER CODE CMD-11 PRICE £105.00 (£126.00)**

## PAG-11c Caller display pager with audible & vibration alert alarm



- Tone/vibration alert carer pager with caller display
- Can pair with up to 99 Rondish transmitters directly
- Sensor type and caller ID assigned during setup
- Powered by 2 x AA batteries (typical 6 months)

**ORDER CODE PAG-11C PRICE £31.00 (£37.20)**

## CMU-02 Wireless nurse call alarm station for up to 999 devices



- Mains powered with battery backup (4 x D-Cell)
- Call and system fault reporting
- Call silence button with delayed alarm repeat
- Displays call type, zone and room/ward ID
- Fully monitored device connection status

**ORDER CODE CMU-02 PRICE £205.00 (£246.00)**



## RONDISH BRAND ADDITIONAL/REPLACEMENT ACCESSORIES

**BMS-03** Non-radio alarm controller for SEN-1/SEN-2/FMAT-2 alarm pads  
**ORDER CODE BMS-03 PRICE £38.00 (£45.60)**



**BMDV-11** Non-radio alarm controller for SEN-1/SEN-2/FMAT-2 alarm pads  
Magnetic pull-away clothing clip alarm and nurse call connection  
**ORDER CODE BMDV-11 PRICE £54.12 (£64.94)**



**DMS-02** Replacement/additional door sensor for DoorWatcher egress alarm  
**ORDER CODE DMS-02 PRICE £248.00 (£297.60)**



**TXP-02R** Replacement/additional resident wristband for DoorWatcher alarm  
**ORDER CODE TXP-02R PRICE £46.00 (£55.20)**



**TM01-V6** Replacement transmitter for cordless bed/chair pads  
**ORDER CODE TM01-V6 PRICE £31.45 (£37.74)**



**WCP-11** Waterproof pull-cord transmitter - use with any Rondish receiver  
**ORDER CODE WCP-11 PRICE £56.00 (£67.20)**



**DRB-11** Wireless call reset for DoorWatcher system  
**ORDER CODE DRB-11 PRICE £68.00 (£81.60)**



**SEN-1** Bed sensor for BMS-03, BMDV11  
**ORDER CODE SEN-1 PRICE £28.00 (£33.60)**



**SEN-2** Chair sensor for BMS-03, BMDV11  
**ORDER CODE SEN-2 PRICE £24.00 (£28.80)**



**FMAT-2** Non-Slip heavy duty large (48" x 24") floor pressure mat  
**ORDER CODE FMAT-02 PRICE £78.00 (£93.60)**



Wire-Free bed and chair occupancy sensor pad with integral transmitter  
Bed pad **ORDER CODE CBP-01 PRICE £70.10 (£84.12)**  
Chair pad **ORDER CODE CCP-01 PRICE £65.80 (£78.96)**



# CALL PENDANT AND SPECIAL SWITCH PRODUCTS

## MPPL-ERTXSET Waterproof call pendant with MPPL pager

Easy to press pendant alarm allows the user to raise an alert to a carer's portable pager, simply by pressing the pendant button.

- Simple wireless SOS call bell solution
- Ideal for the elderly or infirm
- Tone and vibration alert pager
- Pager can pair with up to 9 transmitters
- Pendant is supplied with lanyard and wrist strap
- Portable, lightweight design
- All batteries included



Pager and transmitter kit **CODE MPPL-ERTXSET PRICE £45.00 (£54.00)**

## NMDRX-ERTX Waterproof wireless call pendant with alarm station

The call pendant button signals the alarm station, which is equipped with bright coloured flashing LED's. Useful for hard of hearing or deaf carers.

- Simple to use and set up call solution
- 5-Channel independently coloured flashing LED's
- Output for under pillow vibrating pad
- Front facing and ceiling projected flashing LED's
- Long battery life with low battery warning
- Audible alarm tone with 80-90 dB selection
- Choose your preferred alarm light colour per channel



Call button with alarm **ORDER CODE NMDRX-ERTX PRICE £51.49 (£61.79)**

## CTMV-ERTX Call pendant with 3-channel recordable voice alarm player

Press the waterproof call button to remotely activate playback of your personalised alarm message.

- Useful to generate warnings for blind people
- Create voice reminders or warnings to aid people with dementia
- Create instructional messages to aid people with learning impairment
- Use up to 3 transmitters, each with different voice message
- Single or looped alarm playback



Remote voice alarm kit **ORDER CODE CTMV-ERTX PRICE £49.83 (£59.80)**

## CALL PENDANT AND SPECIAL SWITCH PRODUCTS

### MPPL-PSKIT Soft touch low tactile pillow switch with MPPL pager

Face/head activated switch attaches to pillow to enable a person with limited mobility to activate a help call to a carer pager.

- Soft and comfortable washable switch cover
- Attaches to pillow with included clips
- Signals carer tone/vibration alarm pager
- Transmitter range up to 100M
- Suitable for stroke victims & disabled people
- All batteries included



Pillow switch kit **ORDER CODE MPPL-PSKIT PRICE £164.50 (£196.80)**

### CTMV-PSKIT Soft touch low tactile pillow switch with voice alarm

Face/head activated switch attaches to pillow to enable a person with limited mobility to activate a help call to a carer voice player.

- Use personalised voice recordings as an alarm
- Record up to 3 separate call messages (additional transmitters required)
- Suitable for child carers or people with mild dementia
- Record an instruction for playback when the button is pressed
- Suitable for stroke victims & disabled people
- All batteries included



Pillow switch kit **ORDER CODE CTMV-PSKIT PRICE £169.32 (£203.18)**

### NMDRX-PSKIT Soft touch low tactile pillow switch with alarm station

Face/head activated switch attaches to pillow to enable a person with limited mobility to activate a help call to a carer.

- 5-Channel alarm receiver with 5 selectable alarm colours
- Suitable for hard of hearing/deaf carers
- Audible/flashing light alarm station
- 30 second or looped audible/flashing light alarm
- Transmitter range up to 100M
- Suitable for stroke victims & disabled people
- All batteries included



Pillow switch kit **ORDER CODE NMDRX-PSKIT PRICE £170.99 (£205.19)**



## CALL PENDANT AND SPECIAL SWITCH PRODUCTS

### MPPL-JSK Large surface area call button switch with carer alarm pager

For people that find it difficult or impossible to operate small call pendant buttons, the jelly switch provides a perfect solution. Available with a choice of wireless carer alarm responders.



- Large surface button can be activated with a minimal amount of pressure
- Suitable for patients with poor dexterity
- Image can be placed under the cover to assist patients with dementia or learning difficulties
- Wireless signal to tone/vibration alert alarm pager
- All batteries included

Jelly switch kit with pager **ORDER CODE MPPL-JSK PRICE £82.00 (£98.40)**

Switch with transmitter **ORDER CODE MPPL-JS PRICE £51.49 (£61.79)**

### CTMV-JSK Large surface area call button switch with voice alarm

This jelly switch kit uses the CTMV recordable 3-Channel voice alarm as the carer alert to provide remote controlled playback of voice alarms.



- Record up to 3 separate call messages (additional transmitters required)
- Suitable for child carers or people with mild dementia
- Record an instruction for playback when the button is pressed, "I need help, call an ambulance" etc.

Jelly switch kit with voice alarm **ORDER CODE CTMV-JSK PRICE £86.82 (£104.18)**

Switch with transmitter **ORDER CODE MPPL-JS PRICE £51.49 (£61.79)**

### NMDRX-JSK Large surface area call button switch with carer alarm station

This jelly switch kit uses the NMDRX call station with audible/flashing light alarm.



- 5-Channel alarm receiver with 5 selectable alarm colours
- Suitable for hard of hearing/deaf carers
- 30 second or looped audible/flashing light alarm
- Transmitter range up to 100M
- Suitable for debilitation conditions: Stroke, MS, Parkinson's etc.

Jelly switch kit with alarm station **ORDER CODE NMDRX-JSK PRICE £88.49 (£106.18)**

Jelly switch with transmitter **ORDER CODE MPPL-JS PRICE £51.49 (£61.79)**

# WRIST WORN FALL DETECTION SENSOR WITH CARER ALARM

## MPPL-FA6P Wrist worn gravity activated fall detection alarm

The MPPL-FA6 is an affordable fall detection alarm with touch activated call button. It can be used with an alarm receiver or integrated with one of the Medpage off-site emergency call telephone diallers such as the AD-01.

- Dual fall detection sensor - inertia & gravity
- Use with wrist strap or lanyard (included)
- Transmits alarm signal to carer pager on fall detection
- Comfortable to wear and lightweight
- Touch sensor call button - no need to press button
- Supplied with MPPL carer tone/vibration alert pager
- Suitable for daily and bedtime use - not waterproof



Fall alarm with pager kit **ORDER CODE MPPLFA6P PRICE £58.45 (£70.14)**

Fall alarm only **ORDER CODE MPPL-FA6 PRICE £27.45 (£32.95)**

## TXP11-PAG11 Gravity activated fall detection alarm with pager

Active emergency button/fall detection sends signal to central monitor or included alarm pager

Includes lanyard with anti-ligature point

Optional watch strap (not included)

Accelerometer fall detection sensor

Reset signal can be triggered using magnetic switch

Battery life in excess of 3 years with normal usage

Changeable battery - uses CR2477

Waterproof to IP67



Fall alarm with PAG-11 kit **ORDER CODE TXP11-PAG11 PRICE £63.28 (£75.93)**

Fall alarm only **ORDER CODE TXP-11 PRICE £32.28 (£38.73)**

## NMDRX-FA6P Wrist worn gravity activated fall detection alarm

This kit uses the NMDRX call station with audible/flashing light alarm and the FA6 fall detection wrist band. Fall detection, or touching the call button activates the wireless alarm receiver.

- 5-Channel alarm receiver with 5 selectable alarm colours (example: fall alarm RED bed alarm BLUE)
- Suitable for hard of hearing/deaf carers
- 30 second or looped audible/flashing light alarm
- NMDRX power: 3 x AA batteries or optional AC adaptor



Fall alarm with NMDRX kit **ORDER CODE NMDRX-FA6P PRICE £64.94 (£77.92)**

Fall alarm only **ORDER CODE MPPL-FA6 PRICE £27.45 (£32.95)**

## WIRELESS DOORBELL BUTTON WITH ALARM CHOICES

### NMDRXDB Wireless doorbell button with flashing light alarm receiver

Simple to install and use flashing light doorbell system, ideal for the deaf and hard of hearing.

- Battery operated portable tabletop alarm
- 5-Channel with selectable flashing alarm colour
- Receiver can pair with up to 5 transmitters
- Output for optional under pillow vibrating pad
- 30 second or looped alarm
- Secure doorbell fixings supplied
- All batteries included



Flashing light doorbell kit **ORDER CODE NMDRXDB PRICE £37.49 (£44.99)**  
Doorbell button only **ORDER CODE NMD-DB1 PRICE £8.28 (£9.95)**

### MPPL-DBK Wireless doorbell button with tone/vibration alarm pager

Simple to install and use radio paging doorbell system, ideal for the deaf and hard of hearing.

- Doorbell transmitter button activates pager alarm
- Transmission range up to 80M
- Pager can pair with up to 9 transmitters
- Sturdy pager belt clip
- You can add an NMDRX flashing light (to identify active doorbell signal) alarm or additional pagers



Doorbell with pager **ORDER CODE MPPL-DBK PRICE £31.00 (£37.20)**  
Doorbell button only **ORDER CODE NMD-DB1 PRICE £8.28 (£9.95)**

### CTMV-NDB Wireless doorbell button with voice alarm receiver

A doorbell solution for people with cognitive impairment. Some people fail to respond to a conventional door chime. Simply record a personalised message for playback when the doorbell button is pressed. For example, "Someone is at the front door".

- Can be used with up to 3 transmitters, each with different messages (requires additional transmitters)
- Single or looped message playback



Doorbell with CTMV receiver **ORDER CODE CTMV-NDB PRICE £35.82 (£42.98)**  
Doorbell button only **ORDER CODE NMD-DB1 PRICE £8.28 (£9.95)**



## CALL PENDANT FOR SINGLE AND MULTIPLE USERS

### PAG-11-ULC11 Waterproof call pendant button with caller display alarm pager

The PAG-11 pager is equipped with tone and vibration alert alarms. Transmitters are assigned to a user ID, displayed on the call display.

- Suitable as a single or multi-user call pendant alarm
- Waterproof transmitter with soft touch call button
- Call review button - alarm cancel button - setup menu
- Low battery warning for call button and pager
- PAG-11 is compatible with all Rondish transmitters



Call pendant with pager **ORDER CODE PAG-11-ULC11 PRICE £54.95 (£65.94)**

Call pendant only **ORDER CODE MPPL-ERTX PRICE £14.00 (£16.80)**

### POC-940KIT Long range waterproof call pendant with digital message pager

This call pendant has a typical transmission range of up to 400M. It is compact, lightweight and can be used with either a lanyard or wrist strap, both supplied. The data message pager has tone and vibration alert.

- Suitable for domestic and professional use
- Unlimited call pendants can be monitored by a pager
- Long battery life (up to 6 months) with low battery warning
- Suitable for larger and stone built buildings
- Stores up to 20 call alerts with date and time



Transmitter and pager kit **ORDER CODE POC-940KIT PRICE £104.00 (£124.80)**

Pendant transmitter only **ORDER CODE POC-940 PRICE £46.00 (£55.20)**

### RONWC2 Waterproof pull cord alarm for disabled toilet with emergency call alarm

This is a complete disabled toilet alarm system for use in public and commercial buildings. The battery powered emergency call alarm is wall mounted and has an audible and flashing light alert when signalled by pulling the cord.

- Wireless pull cord transmitter (Range 60M O/F)
- 3M cord length with red pull triangles
- On-device SOS call & wireless reset button
- Install for CQC and DDA compliance



Complete system **ORDER CODE RONWC2 PRICE £168.00 (£201.60)**

Additional pull cord **ORDER CODE WCP11 PRICE £56.00 (£67.20)**

## WIRELESS DOOR ALARMS FOR WANDER DETECTION - SECURITY

### NMDRX-DCTK Medpage wireless door opening alarm with carer alarm receiver

Simple to install door security solution. When the monitored door or window is opened, an alarm signal is transmitted to the portable tabletop alarm receiver, which activates the coloured flashing lights and beeper alarm.

- Dementia and elderly wander detection
- General door and window security
- On device call/panic button
- Sound and visual alert for activated sensor transmitters
- Selectable LED colours to identify active transmitter device
- Powered by 3 x AA batteries



Door transmitter with NMDRX receiver **ORDER CODE NMDRX-DCTK PRICE £49.15 (£58.98)**

### MPPL-DCT Medpage wireless door opening alarm with carer alarm pager

Where a mobile alert to an opening door is required, this kit provides a perfect solution. When the door/window is opened an alarm signal is transmitted to the alarm pager.

- Dementia and elderly wander detection
- General door and window security
- Alarm pager powered by 2 x AA batteries
- Tone/vibration alert with volume selection
- Door alarm powered by 1 x CR2032 battery
- Typical transmitter range 60M



Door transmitter with pager **ORDER CODE MPPL-DCT PRICE £42.66 (£51.19)**

### CTMV-MED-DCT Medpage wireless door opening alarm with voice alarm player

This kit provides a convenient way to remotely trigger playback of a personalised message when a door is opened (internal or external).

- Record safety messages, instructions or warnings for playback when a door is opened
- Useful for care of people with dementia or learning impairment
- Remind people to close a fridge or freezer door
- Useful for blind people to identify opening doors or cupboards



Door transmitter with CTMV receiver **ORDER CODE CTMV-MED-DCT PRICE £47.99 (£57.58)**

# WIRELESS DOOR ALARMS FOR WANDER DETECTION - SECURITY

## PAG-11C-DCT Wireless door alarm with caller display tone/vibrate pager

This kit provides a simple door alarm solution for single or multiple doors. The caller display will show a user ID to indicate door location.

- Suitable for domestic, Care Home and Hospital use
- Use to detect wandering people exiting internal or external doors
- Can integrate with all Rondish wireless care aids
- Pager powered by 2 x AA batteries
- Up to 99 door alarms (or combined devices) can be paired with the pager
- Typical transmitter range 60M. See signal boosters.



Door kit **ORDER CODE PAG-11C-DCT PRICE £42.66 (£51.19)**

## MPPL-DCKITS Wireless door/window security alarm with tone/vibrate pager alarm

Quick and easy to install door/window alarm comprising of MPPL-DC door/window transmitter with our MPPL alarm pager. Simply fix transmitter to door frame using supplied screw pack or Velcro dots. Alarm transmission to pager occurs when door is opened.

- Useful for wander detection
- Door alarm for deaf people
- Vibration and tone alert pager
- Long battery life (typically 6-months)
- Transmitter 1 x 9V battery
- Pager powered by 2 x AA batteries



Door alarm kit **ORDER CODE MPPL-DCKITS PRICE £61.00 (£73.20)**

## POC-DCK Long range (400M) door alarm transmitter with digital message pager

When a door or window alarm transmitter is required with longer transmission, this kit provides a good solution. The compact transmitter (POC-840DC) is programmed to send location messages to the EM300 message displaying pager. Unlimited transmitters can be used with a single or multiple pagers.

- Use for wander detection & general door/window security with door location message
- Suitable for single or multiple alarm transmitter systems
- Suitable for domestic, Care Home and Hospital use
- Factory programmed default message 'Door Alarm 1'
- Tone/vibrate alert pager powered by 1 x AAA battery



Long range door alarm kit **ORDER CODE POC-DCK PRICE £108.00 (£129.60)**



## MEDWFD Wi-Fi connected door alarm with monitoring APP

The MEDWFD is a useful door/window alarm sensor, which connects through your home or business Broadband (2.4GHz only) and is configured via a Smartphone APP (Android or iOS). The APP provides a comprehensive selection of features, including choice of alarm notifications and sensor labelling.



- APP records date & time of door events
- Can also be used on food cupboards and fridges
- Beneficial in the care of people with dementia or cognitive impairment for wander detection
- Can be used with multiple sensors for home and business protection

Wi-Fi door/window sensor **ORDER CODE MEDWFD PRICE £24.95 (£29.95)**

## DWSYS Wristband activated door egress alarm (DoorWatcher)

The DoorWatcher system is designed to allow a person prone to wandering more freedom and independence while in care. The system will alert care staff if the resident approaches a monitored door, (internal or external) corridor, or hazardous area.

The system comprises of a strip style doorway detector that responds to low level signals continually transmitted by the resident wristband. The door strip monitor includes a built in sound alarm with volume control and a transmitter that will signal the nurse call alarm station, CMU-02, CMD-11C or PAG-11C caller display pager.



Installation is very quick and easy requiring no wiring or professional help. Self adhesive pads secure the door strip monitor to the door frame or adjacent wall, the connected power adaptor(s) is then plugged into a mains supply. The resident wristbands are lightweight, comfortable to wear and are powered by sealed in lithium batteries, providing a continuous service life of 36 months.

The wristbands have secure power on/off. Unlimited numbers of wristbands can be used with single or multiple door strip monitors to allow all exit doors to be protected. The advantage of this type of system is the ability to monitor only residents at risk. For a single doorway, one DoorWatcher strip is usually adequate. For a double doorway it is recommended that two DoorWatcher strips are used. The kit price includes all required items for a single resident monitoring system for a single door width.

DoorWatcher kit **ORDER CODE DWSYS PRICE £330.00 (£396.00)**

## MOVEMENT SENSING ALARMS - BED LEAVING - WANDERING

### MPPL-MSKIT PIR movement sensor transmitter with pager

This easy to install movement alarm kit can be used for general security. It is also beneficial to aid care for the elderly and those living with dementia. When the sensor detects movement an alarm signal is transmitted to the pager.

- Use for bed leaving detection and wander alerts
- Protect stairwells and doorways
- Long battery life for PIR & pager (6 months)
- Transmitter range up to 100M
- External power on/off switch for PIR
- Tone/vibration alert pager



PIR with pager **ORDER CODE MPPL-MSKIT PRICE £59.00 (£70.80)**

### NMDRX-PIR2 Movement detection PIR wander detection kit with portable alarm

Detect an elderly patient or a person with dementia wandering during day or night with the Medpage NMDRX-PIR2 wander alarm kit.

- An alert is transmitted to a carer when a person walks through a doorway, enters a room, exits the home or enters a high-risk area
- Uses passive infrared technology to detect movement. The NMDRX is a portable ultra-bright flashing coloured LED alarm receiver with chime alert
- Quick and easy installation
- The PIR sensor can be placed anywhere movement detection is required



PIR with NMDRX receiver **ORDER CODE NMDRX-PIR2 PRICE £51.65 (£61.98)**

### CTMV-MEDPIR2 Movement detection PIR wireless activated voice memo player

This kit has many uses for senior and dementia care. Personal warning messages are recorded to the CTMV alarm receiver. Automatic playback of the message starts when signalled by the wireless motion sensor.

- Record warning or instruction messages for alarm playback
- CTMV has 3-Channels, each with independent voice recordings
- Use additional PIR's for different movement triggered alarms
- PIR re-trigger time settings to prevent constant repeats
- Wireless range up to 60M (see signal repeater)



PIR with CTMV receiver **ORDER CODE CTMV-MEDPIR2 PRICE £49.99 (£59.99)**

## MOVEMENT SENSING BREAK BEAM ALARMS

### FT47C-DTXMK Photoelectric break beam sensor with wireless alarm responder

The FT47C-DTXMK is a compact photoelectric beam set suitable for indoor and outdoor use (shower proof only) to create an invisible narrow beam security solution. An object, person or animal walking through and breaking the beam results in an alarm transmission to the NMDRX alarm responder

- DC powered beam transmitter with battery powered beam receiver reduces cabling and installation time
- Included alarm responder produces a visual and tone notification
- Ideal for short range (200M) security applications, home wander detection.



**ORDER CODE**  
**FT47C-DTXMK**

**PRICE £102.00 (£122.40)**

### FT47C-DTXK Photoelectric break beam sensor alarm kit with long range transmitter and digital message pager

The FT47C-DTXK is a compact photoelectric beam set suitable for indoor and outdoor use (shower proof only) to create an invisible narrow beam security solution. An object, person or animal walking through and breaking the beam results in an alarm transmission to the EM300 data message pager.

- Long range alarm transmission (up to 400M)
- Beam range (strong sunlight) 8-12M (low sunlight) 15-17M
- Programmable alert messages (when used with data pager EM300)
- Working voltage IR transmitter 2 x AA Alkaline batteries
- Typical battery life 12 months 10 ops per day
- Working voltage IR receiver: 12V DC power adaptor or battery
- IP rating: IP55 outdoor use (shower proof)



**ORDER CODE**  
**FT47C-DTXK**

**PRICE £146.00 (£175.20)**



# MOVEMENT SENSING ALARMS - BED LEAVING - WANDERING

## MEDWFM Wi-Fi connected movement sensor with Smartphone alarm APP

Setting up the MEDWFM takes only a minute or two. Download the user APP to your Smartphone and connect to a home or business Broadband router (2.4GHz only). The user-friendly intuitive APP notifies you if a person walks into the detection beam of the PIR, through the internet to your Smartphone, recording the date and time of the alert. You can check if someone is moving around in the home (parents, elderly relative or your own).



- Provides a simple security solution with event date and time record
- Receive an alert if a person enters your home or business premises
- Free to use APP for Android and iOS (no subscription or rental charges)

WiFi movement sensor **ORDER CODE MEDWFM PRICE £29.12 (£34.94)**

## CTM-PIR Battery powered PIR sensor with relay output

The CTM-PIR is a useful battery powered movement detection sensor with relay output. It can be connected to all Medpage/Easylink transmitters and directly to nurse call systems. Use with CTM-3 Voice alarm transmitter for recordable voice prompts.



- Powered by 1 x 9V battery
- Wall fixing swivel bracket

PIR sensor only **ORDER CODE CTM-PIR PRICE £24.25 (£29.10)**

## NMDTX-CTPIRK PIR sensor with POCSAG long-range transmitter & data pager

This high performance kit includes the CTM-PIR coupled to a message programmable POCSAG transmitter (NMDTX) and data message pager (EM300). Movement alarm detections are transmitted to the data message pager (EM300). Suitable for larger buildings and multiple care/security sensor systems.

- Message programmable transmitter with cable connected PIR
- Multi-function transmitter with configurable input (N/O or N/C)
- Transmission range up to 400M



Complete kit **ORDER CODE NMDTX-CTPIRK PRICE £136.25 (£163.50)**

# FLOOR PRESSURE MAT ALARMS

## MPPL-PMK Floor pressure mat alarm with radio pager



Our best-selling floor pressure mat alarm is ideally suited to the domestic user. The mat is easily concealed under a rug or carpet next to an exit door or the user's bed. When the user steps on the mat, the carer is notified via the MPPL pager.



- 100M transmission range
- 3 alert options: vibrate only, vibrate and low tone or vibrate and high tone
- Robust, lightweight pager with a sturdy belt clip
- Low battery warning
- Pressure mat measures 750mm x 450mm with a 3M connection lead

Complete kit **ORDER CODE MPPL-PMK PRICE £75.00 (£90.00)**

## FMAT-2K Large heavy duty non-slip floor pressure mat alarm kit

This kit is more suited to Care Home or Hospital use. The pressure mat is positioned on the floor next to the patient's bed. Standing on the mat results in alarm transmission to the PAG-11C caller display pager. Can be linked into CMD-11 and CMU-02 nurse call stations.



- Mat measures 48" (1016mm) x 24" (609mm)
- Dual surface - carpet grip and rubberised surface
- Includes nurse call transmitter (BTX-04M)
- Includes carer caller display pager (PAG-11C)
- Suitable for single or multi-patient monitoring

Complete kit **ORDER CODE FMAT-2K PRICE £151.00 (£181.20)**

## PM-POCKIT Floor pressure mat alarm with long range (400M) message transmitter & data message pager

The PM-POCKIT is recommended for larger buildings, stone-built properties or buildings with known signalling difficulties. The pressure mat (750mm x 450mm) can be positioned under rugs or carpets and will transmit an alarm signal when a person stands on it.



The kit comprises of NMDTX long range POCSAG transmitter, 1 x CT-3 pressure mat and 1 x EM300 pager. An unlimited number of long range devices can be configured into a multiple care system. Single kit default message "Pressure Mat Alarm".

Complete single kit **ORDER CODE PM-POCKIT PRICE £130.17 (£156.20)**

## FLOOR PRESSURE MAT ALARMS

### WT-ELB Wi-Fi connected floor pressure mat alarm with Smartphone alerts

For families concerned about elderly relatives living alone, the WT-ELB floor pressure mat alarm will provide peace of mind and valuable time/date stamped information of daily living habits. The powerful free Smartphone APP allows you to set times when the sensor is active and choose the type of alert you want to receive to your phone. When activated, the Wi-Fi controlled sensor will send an alarm signal to your Smartphone via your home/business Wi-Fi connection.



- Wander detection - fall prevention
- Controller powered by 2 x AAA batteries
- Under carpet/rug pressure sensor (750mm x 450mm)
- Pair with home/business broadband (2.4GHz only)

Complete kit **ORDER CODE WT-ELB PRICE £86.66 (£103.99)**

### NMDTX-FMATKIT Heavy duty non-slip floor pressure mat with long range transmitter (400M) & digital data message pager

This kit is recommended for larger buildings, stone-built properties or buildings with known signalling difficulties. The pressure mat 48" (1016mm) x 24" (609mm) can be positioned on top of rugs, carpets, or directly on hard floors and will transmit an alarm signal when a person stands on it.



The kit comprises of NMDTX long range POCSAG transmitter, 1 x FMAT-2 pressure mat and 1 x EM300 pager. An unlimited number of long range devices can be configured into a multiple care system. Single kit default message "Pressure Mat Alarm".

Single system kit **ORDER CODE NMDTX-FMATKIT PRICE £183.00 (£219.60)**

### NMDRX-CTMFP Floor pressure mat alarm with transmitter & portable alarm

This kit includes the CT-3 pressure sensor, NMDRX audible/flashing light alarm and the CTM-3 voice alarm transmitter. A warning message can be recorded for playback when a person stands on the mat. A signal is transmitted to activate the NMDRX.



- Dual purpose alarm - voice prompt and carer alert
- Use for dementia care wander detection
- NMDRX can pair with up to 5 sensors to create expanded care systems

Complete kit **ORDER CODE NMDRX-CTMFP PRICE £102.99 (£123.59)**



# BED LEAVING ALARMS FOR SINGLE OR MULTIPLE PATIENTS

## CTM-3BKIT Bed leaving detection alarm with recordable voice alarm and pager

The majority of falls occur in the home when a person leaves their bed. With the CTM-3, a voice recording can be made for playback when a person rises from their bed, along with an alarm transmission to the included MPPL tone/vibrate carer alarm pager.

- Can help reduce falls and potential wandering with advanced warning of bed leaving
- Integral transmitter with operating range of 100M
- Under mattress pressure sensor 750mm x 450mm
- Instant (3 second) or 15 minute delayed alarm
- CTM-3 can switch to floor pressure mat mode
- CTM-3 Powered by 3 x AAA batteries or AC adaptor



Bed leaving kit **ORDER CODE CTM3-BKIT PRICE £94.00 (£112.80)**

## NMDRX-CTM3B Bed leaving detection alarm with wireless audible/flashing alarm

This kit includes the CTM-3 voice alarm transmitter, bed pressure mat CT-3 and the NMDRX wireless tabletop alarm. The NMDRX is a 5-Channel alarm, which enables the inclusion of additional sensors, such as: Door alarms, worn fall sensors (FA6) or PIR (MED-PIR2). A comprehensive home care system is easily configured.

- NMDRX suitable for deaf/hearing impaired carers
- 5 Independent flashing LED colours to assign to alarm channel



Complete kit **ORDER CODE NMDRX-CTM3B PRICE £102.00 (£122.40)**

## NMDTXRXM-CT3B Bed leaving detection kit with channel selectable alarm transmitter

This kit includes the CT-3 bed pressure mat, NMDRX 5-Channel alarm receiver and the NMDTXM multi-function sensor transmitter. The transmitter channel codes are easily changed to allow additional bed monitors to be included in a system (up to 5), with the NMDRX indicating an active alarm with different coloured flashing channel LED's.

- NMDTXM powered by 2 x AAA batteries
- Transmission range up to 200M



Complete kit **ORDER CODE NMDTXRXM-CT3B PRICE £96.31 (£115.57)**

## BED LEAVING ALARMS FOR SINGLE OR MULTIPLE PATIENTS

### NMDTXM-BK Bed leaving alarm kit with under-mattress pressure mat and alarm pager

Place the sensor under the bed mattress, plug the sensor into the controller, switch the controller on. Switch the alarm pager on and you're all set to detect a person leaving their bed.



- Reliable, robust, simple, affordable, and discrete
- All battery operated - typical battery life 12 months
- Alarm transmission range typically 100 Metres
- Selectable alarm delay 3 seconds - 15 minutes
- Mode change from bed to floor pressure mat
- Carer alarm notification pager with Tone/Vibration alert
- Fidget-proof alarm control virtually eliminates false alarms



Complete bed kit **ORDER CODE NMDTXM-BK PRICE £85.32 (£102.38)**

### CBM-03PAGBKIT Bed leaving alarm kit with wire-free under-mattress pressure sensor and caller ID display pager

Totally cordless bed leaving alarm system. The alarm controller transmitter is concealed in the corner of the sensor. With no visible signs of the sensor under the mattress, tampering and damage by wire pulling is reduced. Kit includes PAG-11C caller ID pager and CBP-01 bed sensor with transmitter.



- Suitable for single or multiple bed monitoring
- Compatible with CMU-02, CMD-11C Nurse Call stations



Kit price **ORDER CODE CBM-03PAGBKIT PRICE £101.10 (£121.32)**

### NMDTX-BKIT Under mattress bed leaving sensor for Care Homes

This kit is recommended for larger buildings, stone-built properties or buildings with known signalling difficulties. The bed pressure mat measures 750mm x 450mm and is positioned under the bed mattress. When the patient leaves their bed a data encoded display message is transmitted to the carer pager "off bed". A return to bed is confirmed by new message "on bed".



The kit comprises of NMDTX long range POCSAG transmitter, 1 x CT-3 pressure mat and 1 x EM300 pager. An unlimited number of long range devices can be configured into a multiple care system. Single kit default message "Bed Alarm".

Single user kit **ORDER CODE NMDTX-BKIT PRICE £140.82 (£168.98)**

# BED LEAVING ALARMS FOR SINGLE OR MULTIPLE PATIENTS

## CTMV-CT3 Bed leaving detection alarm with recordable voice prompt

This kit is suited to people who may be confused when waking during the night. Upon rising from the bed a personalised message is played back as a warning or instruction. The CTMV is a receiver with 3 voice alarm channels that can be triggered by remote transmitters. A door alarm, for example, could remind a person a door is open before getting into bed.



- Under mattress pressure sensor 750mm x 450mm
- 3-Channel voice recordings of up to 15 seconds each
- Voice alarm CTMV powered by 3 x AAA batteries (included)

Complete kit **ORDER CODE CTMV-CT3 PRICE £56.83 (£68.20)**

## BMDV-11BK Dual alarm bed monitor with recordable voice prompt and clothing clip pull-away alarm

This multi-function bed/chair exit alarm monitor from Rondish, offers state-of-the-art fall prevention technology for home and professional care use. The dual alarm trigger function (sensor pad or clothing clip) provides options for bed, chair or W/C. Can be cable connected to Nurse Call system.



- Pull cord clip alarm - use for bed, chair, W/C
- Bed under-mattress pressure sensor
- Recordable voice prompt - "stay in bed" etc.

Complete kit **ORDER CODE BMDV-11BK PRICE £82.12 (£98.54)**

**Rondish**

## BMS-03BK Bed leaving alarm under-mattress pressure sensor, audible alarm and nurse call connection output

Basic bed leaving alarm which produces an audible alarm when a patient rises from their bed.



- Under mattress pressure sensor
- Sensor controller with audible alarm
- Sensor pad fault/unplugged alarm
- Powered by 3 x AAA batteries
- Nurse Call connection output jack socket

Controller with sensor **ORDER CODE BMS-03BK PRICE £66.00 (£79.20)**



## CHAIR LEAVING ALARMS FOR SINGLE OR MULTIPLE PATIENTS

### CTM3-CKIT Chair leaving detection alarm with recordable voice alarm and pager

Falls often occur in the home when a person rises from their chair. With the CTM-3, a voice recording can be made for playback when a person rises from their chair, along with an alarm transmission to the included MPPL tone/vibrate carer alarm pager.

- Can help reduce falls and potential wandering with advanced warning of chair leaving
- Integral transmitter with typical range of to 100M
- Under cushion pressure sensor 380mm x 300mm
- Powered by 3 x AAA batteries or optional AC adaptor
- Monitor holding bracket included



Complete kit **ORDER CODE CTM3-CKIT PRICE £90.75 (£108.90)**

### NMDRX-CTMCK Chair leaving detection alarm with wireless audible/flashing alarm

This kit includes the CTM-3 recordable voice alarm transmitter, chair pressure mat CT-2 and the NMDRX wireless tabletop alarm. The NMDRX is a 5-Channel alarm, which enables the inclusion of additional sensors, such as: Door alarms, worn fall sensors (FA6) or PIR (MED-PIR2). A comprehensive home care system is easily configured (See page 12).

- NMDRX suitable for deaf/hearing impaired carers
- 5 Independent flashing LED colours to assign to alarm channel - 100M wireless range
- Voice reminder and alarm receiver alert on chair exit



Complete kit **ORDER CODE NMDRX-CTMCK PRICE £99.99 (£119.99)**

### NMDRX-NMDTXCKIT Chair leaving detection kit with channel selectable alarm transmitter

This kit includes the CT-2 chair pressure mat, NMDRX 5-Channel alarm receiver and the NMDTXM multi-function sensor transmitter. The transmitter channel codes are easily changed to allow additional chair monitors to be included in a system (up to 5) with the NMDRX indicating an active alarm with different coloured flashing channel LED's.

- NMDTXM powered by 2 x AAA batteries
- Transmission range up to 100M
- Advanced warning fall prevention/wander detection
- Portable alarm receiver with audible/flashing LED's



Complete kit **ORDER CODE NMDRX-NMDTXCKIT PRICE £88.81 (£106.57)**

## CHAIR LEAVING ALARMS FOR SINGLE OR MULTIPLE PATIENTS

### NMDTXM-CK Chair leaving alarm kit with under-cushion pressure mat and alarm pager

Place the sensor under the chair cushion, plug the sensor into the controller, switch the controller on. Switch the alarm pager on and you're all set to detect a person leaving their chair.

- Reliable, robust, simple, affordable, and discrete
- All battery operated - typical battery life 12 months
- Alarm transmission range typically 100 Metres
- Selectable alarm delay 3 seconds - 15 minutes
- Mode change from chair to floor pressure mat
- Carer alarm notification pager with Tone/Vibration alert
- Fidget-proof alarm control virtually eliminates false alarms

Complete chair kit **ORDER CODE NMDTXM-CK PRICE £82.32 (£98.78)**



### CBM-03PAGCKIT Chair leaving alarm kit with wire-free under-cushion pressure sensor and caller ID display pager

Totally cordless chair leaving alarm system. The alarm controller transmitter is concealed in the corner of the sensor. With no visible signs of the sensor under the cushion, tampering and damage by wire pulling is reduced. Kit includes PAG-11C caller ID pager and CCP-01 chair sensor with transmitter.

- Expandable from a single to multiple patient system
- Alarm signal transmitted to pager upon rise from chair
- Compatible with CMU-02 & CMD-11C nurse call stations

Kit price **ORDER CODE CBM-03PAGCKIT PRICE £91.80 (£110.16)**



### NMDTX-CK Under-cushion chair leaving alarm for Care Homes

This kit is recommended for larger buildings, stone-built properties or buildings with known signalling difficulties. The chair pressure mat measures 380mm x 300mm and is positioned under the chair cushion. When the patient rises from their chair a data encoded display message is transmitted to the carer pager. For example, 'Chair Alarm 1'.

The kit comprises of NMDTX long range POCSAG transmitter, 1 x CT-2 pressure mat and 1 x EM300 pager. An unlimited number of long range devices can be configured into a multiple care system. Single kit default message "Chair Alarm".

Single user kit **ORDER CODE NMDTX-CK PRICE £127.00 (£152.40)**



# CHAIR LEAVING ALARMS FOR SINGLE OR MULTIPLE PATIENTS

## CTMV-CT2 Chair leaving detection alarm with recordable voice prompt

This kit is suited to people who may be unsteady on their feet when rising from their chair. A personalised recordable message is played back as a warning or instruction.



- Under cushion pressure sensor 380mm x 300mm
- 3-Channel voice recordings of up to 15 seconds each
- Integral alarm receiver with Nurse Call connection output
- Voice alarm CTMV powered by 3 x AAA batteries (included)
- Advanced warning of potential fall or wandering
- Single or looped message playback

Complete kit **ORDER CODE CTMV-CT2 PRICE £56.83 (£68.20)**

## BMDV-11CK Dual alarm chair monitor with recordable voice prompt and clothing clip pull-away alarm

This multi-function chair exit alarm monitor from Rondish, offers state-of-the-art fall prevention technology for home and professional care use. The dual alarm trigger function monitor provides options for under chair cushion (pad) or pull-off clothing clip, also useful for W/C. Can be cable connected to Nurse Call.



- Pull cord alarm attaches to clothing
- Under chair cushion pressure sensor
- Recordable voice prompt - "stay in your chair" etc.

Complete kit **ORDER CODE BMDV-11CK PRICE £78.40 (£94.08)**

## BMS-03CK Chair leaving alarm under-cushion pressure sensor, audible alarm and nurse call connection output

Basic chair leaving alarm which produces an audible alarm when a patient rises from their chair.



- Under-cushion pressure sensor
- Sensor controller with audible alarm
- Volume adjustment
- Sensor pad fault/unplugged alarm
- Powered by 3 x AAA batteries
- Nurse Call connection output jack socket

Controller with sensor **ORDER CODE BMS-03CK PRICE £62.00 (£74.40)**



## OFF-SITE ALERT HOME CARE SYSTEMS - ACTIVITY/INACTIVITY

### **CTMVADFA6 Wrist worn fall detection alarm with touch sensor panic button, CTMV voice alarm receiver and telephone emergency dialler**

This system helps to ensure the safety of a person living independently at home and provides peace of mind to distant carers. The wrist worn fall alarm is configured to detect drop falls, then transmits an alarm signal to the CTMV voice alarm, which triggers dialling of your stored emergency contacts.



The CTMV will announce your recorded message repeatedly until reset, enabling call recipients to hear the message detected by the dialler intercom microphone. Requires telephone land line.

Complete system **ORDER CODE CTMVADFA6 PRICE £109.11 (£130.93)**

### **BC346-EM300 Waterproof fall detection bracelet with long range call button**

The waterproof BC346 has integral gravity and inertia sensors for fall detection and a convenient easy-press SOS call button. This kit includes a digital message pager with tone and vibration for alarm notification alerts. Ideal for Senior, Epilepsy and Parkinson's care.

The transmitter has a range of up to 400M making it suitable for domestic and Care Home use. The optional Wi-Fi connected control hub (MP-900) provides off-site alarm alerts to Smartphone and setup of any sensors added to the kit (P.11) including: door alarms, bed & chair occupancy sensors. Kit includes: 1 x fall sensor BC346, 1 x EM300 pager and all required batteries.



**BC346-EM300**

**ORDER CODE BC346-EM300 PRICE £139.00 (£166.80)**

### **ELK4-PO GSM (mobile network) connected SOS help call dialler**

The ELK4-PO uses the mobile telephone network to automatically call emergency contacts when the pendant button is pressed. The system can be used at home and as it uses the mobile network it can be taken anywhere, on holiday etc.

- Mobile network connected emergency call pendant alarm kit
- Large SOS button & included portable call pendant initiates emergency contact dialling
- 2-Way voice intercom - Room temperature monitoring
- Intuitive Smartphone APP allows quick & easy setup of all functions
- No monitoring or subscription charges
- Requires mobile SIM card (not included)



Complete kit including batteries **ORDER CODE ELK4-PO PRICE £49.95 (£59.94)**